


HOW THE VERIFLEX™ PIOL WORKS


To understand how the VeriFlex™ pIOL works, you first have to understand how the eye works. There are structures in your eye needed to convert light into the images you see. When light enters the eye, your cornea and crystalline lens focus it onto the retina, which transforms it into the electrical signals that are sent to your brain. Because the cornea has the greatest focusing power, it's the cornea's shape that determines the majority of the eye's focusing ability.

BEFORE VeriFlex™ pIOL



In normal eyes, the cornea bends light rays to focus directly on your retina. In nearsighted eyes, the cornea's shape is too steep (or the eye is too long), and light is focused in front of your retina. This causes blurry vision.

AFTER VeriFlex™ pIOL



The VeriFlex™ pIOL is designed to permanently treat your nearsightedness by correcting the way light focuses on the retina. The VeriFlex™ pIOL sits behind your cornea, and redirects light so you can enjoy clearer vision.

NOW IT'S YOUR TURN FOR PRECISION VISION

LOST WITHOUT YOUR GLASSES? VERIFLEX™ PHAKIC IOL CAN GIVE YOU A CLEARER VIEW OF LIFE.

For very nearsighted people, procedures like Custom LASIK may not be an option, leaving you stuck with thick glasses or contact lenses — until now. The VeriFlex™ phakic intraocular lens (pIOL) provides “precision vision,” a high-tech, yet simple procedure that not only can improve your vision, but also your life.

What is the VeriFlex™ pIOL?

The VeriFlex™ pIOL has been demonstrated to correct moderate to severe nearsightedness. The procedure involves placing the VeriFlex™ pIOL inside your eye to provide high-quality, high-definition vision without the need for glasses or contact lenses in most situations. Although the VeriFlex™ pIOL is intended to be permanent, the lens can be removed if desired.

* Market Scope 2005 Worldwide IOL Market Report.

The VeriFlex™ Phakic IOL is indicated for the correction of axial myopia, (range -2.0D to -14.5D in 0.5 diopter increments) for patients 18 to 60 years of age, when there are no compromising ocular pathology(ies) and the patients are in otherwise good general health. For a complete list of warning, precautions and adverse events, please refer to the package insert. Rx Only.

VeriFlex and the VeriFlex logo are trademarks of Advanced Medical Optics, Inc.

April 2006



Very precise vision for nearsighted people who want independence from glasses



IS THE VERIFLEX™ PIOL RIGHT FOR YOU?

If you are very nearsighted and tired of contact lenses and thick glasses, and at least 18 years of age, you may be a candidate for the VeriFlex™ phakic IOL. Your doctor will perform a series of tests and discuss options with you to help you decide if you meet the following criteria for having a VeriFlex™ pIOL procedure:

- Your eyes are healthy with no eye disease; iris, pupil or corneal abnormality; or history of retinal detachment.
- Your nearsightedness can be corrected with lenses between -2.0 to -14.5 diopters.
- You have 2.5 diopters or less of astigmatism.
- Your eyes are stable and your nearsightedness has not changed more than 0.50 diopters for 1 year prior to surgery.
- You are not currently pregnant or nursing.

WHAT IS A “PHAKIC INTRAOCULAR LENS”?

The VeriFlex™ phakic IOL is a lens placed behind your cornea. The word “phakic” means that your natural crystalline lens is left in the eye (not removed like in cataract surgery). This is important because your natural lens helps your eye adjust between seeing objects that are near and far.

A PROVEN HISTORY AND HIGH SATISFACTION

The VeriFlex™ pIOL is part of a family of iris-attached lenses used successfully in more than 160,000 procedures worldwide.* The VeriFlex™ pIOL design provides precise and predictable correction for moderate to severe myopia.


ABOUT THE VERIFLEX™ PIOL PROCEDURE

The VeriFlex™ phakic IOL implantation is an outpatient treatment that usually takes about 15 to 30 minutes.

One eye is treated at a time. Generally, the procedure goes as follows:

- 1 Eye drops will be given to reduce your pupil size.
- 2 For better access to your eye, your doctor will use an instrument to hold your eyelids open during the procedure.
- 3 A local and/or IV anesthetic is given, so the procedure is virtually painless.
- 4 A small incision is made in the cornea for the VeriFlex™ pIOL to be inserted in the space between the iris and the cornea.
- 5 The VeriFlex™ pIOL is centered in front of the pupil (the black part in the center of your eye), and is gently attached to the iris (the colored part of your eye) to hold the lens in place.
- 6 A temporary shield will be placed over your eye to protect it during the time just after the procedure.

HOW THE VERIFLEX™ PHAKIC IOL IS PLACED IN YOUR EYE



Through a small incision, your doctor places the VeriFlex™ phakic IOL behind your cornea. You generally do not feel or see the lens once it is in place.

Only your doctor can determine if you are a candidate for the VeriFlex™ pIOL. Be sure to review the complete Patient Information Booklet (available through your doctor) for more details, and discuss its contents with your doctor so that you have all your questions answered.

